

Product Catalogue

2009-06-01



The Art of Mobility™

www.permobil.com

Contents

C300	3
Accessories, C300 chassis	4
C400	5
Accessories, C400 chassis	6
C500	7
Accessories, C500 chassis	8
C350	9
Accessories, C350 chassis	10
Street.....	11
X850	12
Accessories, X850 chassis	13
K300	13
Chairman HD	14
Accessories, Chairman HD.....	15
Conversion articles, Chairman HD	16
K300	17
Accessories, K300 chassis	18
K450	19
Accessories, K450 chassis	20
Corpus seat.....	21
Conversion articles, Corpus	24
Seat and back cushions, Corpus	25
VS seat	26
Conversion articles, VS.....	28
Seat and back cushions, VS	29
PS seat.....	30
TS seat	31
Conversion articles, TS.....	32
CS seat	33
Miniflex	34
Conversion articles, Miniflex	35
UniTrack.....	36
Support, UniTrack	38
Spare parts, UniTrack	40
Control system.....	42
DriveLink MultiController	44
DL Accessories	46
Permolock C	47
Permolock C3	48

Permobil AB
Box 120, 861 23 Timrå, Sweden
+46-60-59 59 00

The latest version of this catalogue can be downloaded from our website as a PDF file. Our website also contains our other products, news and other information.

www.permobil.com



C300

A flexible, small combination chassis

The most economical solution

Strong outdoors despite its small size

Four seat options

Seat tilt Seat lift Leg rest Back rest



Turning radius
625 mm



Seat height
510-710 mm



Turning circle
1030 mm



C300 Corpus

Ergonomically approved

Modules allow for changes to be made

Can be fitted with all accessories and
additional functions

Rnet Electronics



Turning radius
625 mm



Seat height
460-660 mm
510-710 mm



C300 TS

Built-in UniTrack

Seamless adjustment of seat depth and
arm rest height

Functional design

VR2 Electronics



Turning radius
625 mm



Seat height
510-710 mm



C300 PS

A smart seat that is easy to adjust

A sliding, telescopic frame and expanding seat
pan

VR2/Rnet Electronics



Turning radius
625 mm



Seat height
490-690 mm



C300 CS

Basic but comfortable seat with manual
adjustment of backrest angle

VR2 Electronics

Max. speed: 7 km/h
 Range: 30-35 km
 Obstacle clearance ability: 70 mm
 Colors available: metallic red, metallic blue, metallic green, metallic black

Accessories, C300 chassis



Battery, 60 Ah
101969-99-0



Charger, JK Medico 10 A
104750-99-0



Lights
103601-99-0, VR2
103564-99-0, Rnet




Bag
104753-99-0



Solid wheel
103499-99-0



Cover
101888-99-0



Permolock C3
1823300



C400

A strong, flexible combination chassis

Drives well indoors and outdoors

An intermediate chassis for slightly smaller areas

Seat tilt Seat lift Leg rest Back rest



Turning radius
665 mm



Seat height
520-720 mm
Lowrider
460-610 mm



Turning circle
1075 mm



Two seat options

C400 Corpus

Ergonomically approved

Modules allow for changes to be made

Can be fitted with all additional functions

Available in a junior size

Rnet Electronics

Standing



Turning radius
665 mm



Seat height
530-780 mm



Turning circle
1025 mm



C400 VS

The standing function is both useful and practical

Can be raised from various positions with one button

Favorite positions are saved to make life easier

Also available in a junior size

Rnet Electronics

Max. speed:	7 km/h
Range:	30-35 km
Obstacle clearance ability:	70 mm
Vehicle lock:	Permolock C
Colors available:	metallic red, metallic blue, metallic green, metallic black

Accessories, C400 chassis



Battery, 73 Ah
101989-99-0



Luggage basket
101500-99-0



Lights
103635-99-0



Seat stabilizer
100570-99-0



Solid wheel
103889-99-0



Permlock Plunge
102842-99-0



Charger, JK Medico 10 A
104750-99-0



Bag
104753-99-0




Running time meter set
101957-99-0



Support wheel set
102210-99-0



Transport brackets
101515-99-0



Permlock C
102840-99-0




C500

A strong combination chassis that drives well outdoors

Available as a Lowrider model for a lower seat height

Or as an S model, a faster version

Seat tilt Seat lift Leg rest Back rest



Two seat options

C500 Corpus

Ergonomically approved

Modules allow for changes to be made

Can be fitted with all additional functions

Available in a junior size

Rnet Electronics

Standing



C500 VS

The standing function is both useful and practical

Can be raised from various positions with one button

Favorite positions are saved to make life easier

Available in a junior size

Rnet Electronics

Max. speed:	7.5 km/h
With S pack:	10 km/h
Range:	30-40 km
Obstacle clearance ability:	70 mm
Vehicle lock:	Permlock C
Colors available:	metallic red, metallic blue, metallic green, metallic black

Accessories, C500 chassis



Battery, 73 Ah
101989-99-0



Charger, JK Medico 10 A
104750-99-0



Luggage basket
101500-99-0



Trelleborg tire
101378 -99-0 (coarser tread)
Solid wheel
103499-99-0



Lights
103635-99-0



Support wheel set
102210-99-0



Seat stabilizer
100570-99-0



Transport brackets
101515-99-0




Cover
101888-99-0



Permlock Plunge
102842-99-0



Running time meter set
101957-99-0



Permlock C
102840-99-0



C350

A solid back wheel drive chassis

Strong and easy manoeuvred

Smart design that shows

Seat tilt Seat lift Leg rest Back rest



Turning radius
890 mm



Seat height
510-710 mm



Three seat option

C350 Corpus

Ergonomically approved

Modules allow for changes to be made

Can be fitted with all accessories and additional functions

Rnet Electronics



Turning radius
890 mm



Seat height
460-660 mm
510-710 mm



C350 TS

Built-in UniTrack

Seamless adjustment of seat depth and arm rest height

Functional design

VR2 Electronics



Turning radius
890 mm



Seat height
510-710 mm



C350 PS

A smart seat that is easy to adjust

A sliding, telescopic frame and expanding seat pan

VR2/Rnet Electronics

Max speed: 10 km/h
 Range: 30-35 km
 Obstacle clearance ability: 70 mm
 Colors available: metallic red, metallic blue,
 metallic yellow, metallic black
 metallic silver

Accessories, C350 chassis




Battery, 60 Ah
101969-99-0



Lights
103035-99-0



Solid wheel
103499-99-0



**Locking pins,
Permolock C3**
1824229



Charger, JK Medico 10 A
104750-99-0



Bag
104753-99-0



Cover
101888-99-0

Permolock C3
1823300



Street

Back wheel drive combination chassis

Fast and flexible in town

Does 12 km/h

Seat tilt Seat lift Leg rest Back rest



Turning radius
900 mm



Seat height
560-810 mm



Turning circle
1220 mm



Street Corpus

Ergonomically approved

Can be fitted with all additional functions

Modules allow for changes to be made

Available in a junior size

Rnet Electronics

Max. speed:

12 km/h

Range:

25-35 km

Obstacle clearance ability:

70 mm

Color options:

metallic brown, gold/black,
metallic blue, silver/black



Battery, 73 Ah

101989-99-0



Charger, JK Medico 10 A

104750-99-0

Bag

104753-99-0



Lights

102058-99-0



X850

A powerful outdoor chassis that drives everywhere

Fast, strong and robust

Decent suspension for tougher surfaces

Seat lift Back rest Seat rotation Adjustable wheelbase Seat tilt



Turning radius
1680 mm



Seat height
560-760 mm



Turning circle
1930 mm



Three seat options

X850 Corpus

Ergonomically approved

Modules allow for changes to be made

Can be fitted with a number of additional functions

Available in a junior size

Safegate Electronics



Turning radius
1680 mm



Seat height
530-730 mm



Turning circle
1930 mm



X850 Miniflex

A child seat for children up to 35 kg

Easy to use so that the child feels safe

Safegate Electronics



Turning radius
1680 mm



Seat height
550-750 mm



Turning circle
1930 mm



X850 CS

Basic but comfortable seat with manual adjustment of backrest angle

Powered or manual seat rotation for easier transfers (option)

Safegate Electronics

Max. speed: 15 km/h
 Range: 35-45 km
 Obstacle clearance ability: 110 mm
 Color options: metallic red, metallic blue, metallic black

Accessories, X850 chassis



Tire, fine tread
103016-99-0

Tire, coarse tread
103017-99-0



Battery, 97 Ah
101999-99-0



Charger, JK Medico 10 A
104750-99-0

Bag
104753-99-0



Luggage basket
103019-99-0



Luggage box
103039-99-0 Front opening
103038-99-0 Rear opening




Running time meter set
101957-99-0



X850 footplate
101892-99-0



Crutch holder
100292-99-0 Right
100293-99-0 Left



Length adjustment, electric
103012-99-0



Cover
101899-99-0 X850 With frame
101888-99-0 X850 Without frame



Frame
103171-99-0 Miniflex
103172-99-0 CS



Electric seat rotation
103009-99-0
Manual seat rotation
103023-99-0



Mirror, incl. bracket
103101-99-0 Right
103137-99-0 Left



Chairman HD

A robust combination chassis for bigger users

Supports weights up to 200 kg

Reinforcement in the seat frame, back and arm rest



Chairman HD Corpus

Ergonomically approved

Larger than a standard Corpus

Reinforced in the frame, back and arm rest

Can be fitted with all additional functions

Accessories, Chairman HD chassis



Battery, 73 Ah
101989-99-0



Charger, JK Medico 10 A
104750-99-0



Running time meter set
101957-99-0



Rocker set, switch panel
100524-99-0



Belt with quick-release lock, HD
100507-99-0



Belt, self-adjusting, long
100600-99-0

Max. speed: 6.5 km/h
 Range: 30-35 km
 Obstacle clearance ability: 60 mm
 Color options: black

Accessories, Chairman HD



Head rest, adjustable

104130-99-0 Fabric
 104131-99-0 Imitation leather

Lateral head rest adjustment

103697-99-0



Calf rest, incl. bracket

100514-99-0 Fabric
 100515-99-0 Imitation leather



Thigh support, complete, long

100512-99-0 Fabric
 100513-99-0 Imitation leather



Thigh support, complete, short

100510-99-0 Fabric
 100511-99-0 Imitation leather



Arm rest bag

100671-99-0 Right
 100672-99-0 Left
 102588-99-0 Front of seat frame



Mug/Phone holder

103044-99-0



Table, angle adjustable

103426-99-0 (500x332)



Co-Pilot

103181-99-0



Parallel movement panel holder, Pilot+

100819-99-0 Right
 100824-99-0 Left



Cover

101888-99-0



Foot strap with plastic lock

103427-99-0



Foot strap with padding

103104-00-0

Conversion articles, Chairman HD




**Seat reduction,
Chairman HD**
309138-28-0




**Seat extension for
585 depth**
309253-28-0



**Seat extension for
635 depth
Chairman HD**
309254-28-0 Seat extension +100



**Seat widening for
535 depth**
309150-28-0



**Seat widening for
585 depth**
309238-28-0



**Seat widening for
635 depth**
309239-28-0



Back widening set
309223-40-0



**Trunk rest, HD,
narrow back**
309217-00-0



**Support for the curve
of the back**
306396-99-0



**Thigh rest bracket/
Table bracket**
301612-28-0



Arm Rest Neutral, Fabric
1822784, 350 mm
1822785, 400 mm



**Arm Rest Neutral,
Imitation Leather**
1822788, 350 mm
1822789, 400 mm



K300

A chair to grow in

Cool kids in all sizes can get in on the action

Seat tilt Seat lift



Turning radius
625 mm



Seat height
420-620



Turning circle
1146 mm



K300 PS Junior

A child seat for children up to 76 kg.
Adjustable seat pan fits kids of all sizes
Rnet Electronics

Max. speed:

8 km/h

Range:

30-35 km

Obstacle clearance ability:

60 mm

Colors available:

chilipepper red, chrystal blue
black onyx, platinum silver,
popstar pink, titan blue



Charger, JK Medico 10 A
101958-99-0



Battery, 60 Ah
101969-99-0



Solid wheel
103499-99-0




Lights
103564-99-0



Permolock C3
1823300



Cover
101888-99-0



**Locking pins,
Permolock C3**
1824229

Max. speed: 7 km/h
 Range:: 30-35 km
 Obstacle clearance ability: 70 mm
 Colors available: metallic red, metallic pink, metallic blue, chrystal blue, metallic black metallic green, metallic silver,

Accessories K300



Joystick Compact Advance
103814-99-0




Back box incl. bag hook set
101888-99-0



Headrest
103395-99-0

Adapter head rest slingback
104148-99-0



Bag hook set
101958-99-0



Table, angle adjustable
103426-99-0



Calf support
102981-99-0, flat
103058-99-0, large
103060-99-0, x-large
102978-99-0, medium
103055-99-0, small
Bracket
103678-99-0



Belt buckle, manual adjustm.
100518-99-0



Trunk support
103056-99-0




Belt, self-adjusting
100519-99-0



Push handle
103811-99-0



Belt, self-adjusting, long
100600-99-0



Quick lock, backrest adjustm.
104616-99-0



K450

Broad and solid tyres for a smooth and comfortable ride
Rear wheel drive and an full (all wheel) suspension system
Unique seat that can take you from an up-right position down to the floor level in a few second

Seat tilt Seat lift Legrest Backrest



Turning radius
1060 mm

Seat height
120-650 mm



K450 MX

The K450 MX is for users weighing up to 60 kg, with a flexible seat width that can be adjusted from 25-40 cm. Perfect for growing children.
Rnet Electronics



Battery, 60 Ah
101989-99-0




Charger, JK Medico 10 A
104750-99-0




Cover
101888-99-0



Lights
104540-99-0



Locking pins, Permlock C3
1824229



Transport bracket
104705-99-0



Permlock C3
1823300

Max. speed: 8,5 km/h
 Range: 25-30 km
 Obstacle clearance ability: 60 mm
 Colors available: metallic red, metallic pink,
 metallic blue, metallic black,
 metallic silver; yellow

Accessories K450



Joystick Compact Advanced
103814-99-0



Footstrap with plastic lock
103427-99-0

Footstrap with padding
103104-99-0



Co-Pilot, attendant control
104571-99-0



Table adjustable angle
100575-99-0



Headrest, 3 models
104572-99-0, Imitation leather
104573-99-0, fabric
104339-99-0, Imitation leather small

Lateral headrest adjustment
104344-99-0, 250 mm
104345-99-0, 300-350 mm



Belt buckle, manual adjustment
100518-99-0



Thigh support
103081-99-0, large
103077-99-0, småö



Trunk support
104574-99-0, medium
104575-99-0, small

Corpus seat

As a very ergonomic solution and with its electric seat angle and back rest angle adjustment, Corpus is the most popular electric seat system in the world. Corpus was designed by one of the leading ergonomic experts in the world to fit the natural contours of the body. The seat's modular design improves functionality by means of changes in size and

electrical functions. The customization potential of the Corpus seat system is virtually unlimited. The depth, width, height and angle of several seat components can be adjusted or adapted. There is a wide range of accessories for the seat, including the UniTrack accessory system.



**Parallel movement
panel holder**
104302-99-0



Table, angle adjustable
103426-99-0 (500x332)



Table, adjustable angle
100575-99-0 (350 mm)
100576-99-0 (400 mm)
100577-99-0 (450 mm)



Small table, incl. bracket
100578-99-0 Fitted on the right
101691-99-0 Fitted on the left
(250x250 mm)



**Table, swiveling with side
cutout**
100571-99-0 (600 mm)
100572-99-0 (650 mm)

With R-net Electronics
103982-99-0
103983-99-0




**Table, swiveling with
central cutout**
100573-99-0 (600 mm)
100574-99-0 (650 mm)

With R-net Electronics
103980-99-0 (600 mm)
103981-99-0 (650 mm)



Back box
102059-99-0



Bag hook set
101473-99-0



Arm rest bag
100671-99-0 Right
100672-99-0 Left
102588-99-0 Front of seat frame



**Belt buckle,
manual adjustment**
100518-99-0



Belt, self-adjusting
100519-99-0



Belt, self-adjusting, long
100600-99-0

Corpus seat



Incontinence protection
102806-99-0 (480 width, 510 and 560 depth)
102807-99-0 (430 width, 460, 510 and 560 depth)



For Comfort seat cushions, see Corpus Seat and Back Cushions



Incontinence protection
102808-99-0 (430 width, 360 and 410 depth)
102809-99-0 (480 width, 460 depth)



Comfort seat plate
101227-99-0



Push handle set
100975-99-0



Arm rest support handle
103194-99-0



Support handle set
103651-99-0



Arm rest raising set
100523-99-0



Side support
101261-99-0 Right outside
101262-99-0 Left outside



Arm rest reinforcement set
101510-99-0



Adjustable head rest
103460-99-0

Lateral headrest adjustment
103697-99-0



Head rest, adjustable
104130-99-0 Fabric
104131-99-0 Imitation leather



Mirror, incl. bracket
101440-99-0



Inflatable support for the curve of the back
101972-99-0



Co-Pilot
103181-99-0

With R-net Electronics
104003-99-0



Oxygen holder
101964-99-0

Corpus seat



Mug/Phone holder
103044-99-0



Crutch holder
101956-99-0



Foot strap with plastic lock
103427-99-0



Foot strap with padding
103104-00-0



Leg divider
100978-99-0

Support for Corpus; see UniTrack accessory system.

Conversion articles, Corpus



Arm Rest Neutral, Fabric
1822782, 250 mm
1822783, 300 mm
1822784, 350 mm
1822785, 400 mm




Arm Rest Neutral, Imitation Leather
1822786, 250 mm
1822787, 300 mm
1822788, 350 mm
1822789, 400 mm




footplate, large
305592-28-0



Divided footplate
306614-28-1 Left
306614-28-2 Right



Seat reduction for seat size D360 W430 without cushion
306740-28-UD-0 Reduction -100



Seat reduction set for size D410 W430 without cushion
305265-28-0 Reduction -50



Seat extension, 560 mm
305260-28-UD-0 Without cushion



Seat extension, 510 mm
305255-28-0 Extension +50



Seat widening, 480 mm
1822383 (460 mm depth)
1822384 (510 mm depth)
1822385 (560 mm depth)



Back set, 535 mm, low
310027-00-0 (340 mm wide)
310030-00-0 (390 mm wide)
310033-00-0 (440 mm wide)



Back set, slide
305200-00-UD-0 (390 mm wide)
305252-00-UD-0 (440 mm wide)



Back set, 605 mm, STD
305480-00-UD-0 (340 mm wide)
305300-00-UD-0 (390 mm wide)
305490-00-UD-0 (440 mm wide)



Support for the curve of the back
306396-99-0



Shaft, PLG arm rest bracket
302722-99-0 for 340 and 390 mm back width
304968-99-0 for 440 mm back width

Seat and back cushions, Corpus

The dimensions below are for the seat. The actual dimensions of the cushion are normally a little bigger.



BACK CUSHIONS CORPUS

Check	Mesh	Imitation leather	Width x height (seat)	Width x height (cushion)
307506-99-80-0	307506-99-60-0	307506-99-90-0	340x535	370x530
307507-99-80-0	307507-99-60-0	307507-99-90-0	390x535	420x530
307508-99-80-0	307508-99-60-0	307508-99-90-0	440x535	470x530
306475-99-80-0	306475-99-60-0	306475-99-90-0	340x605	370x600
306430-99-80-0	306430-99-60-0	306430-99-90-0	390x605	420x600
306473-99-80-0	306473-99-60-0	306473-99-90-0	440x605	470x600

SEAT CUSHIONS CORPUS

Check	Mesh	Imitation leather	Width x height (sits)	Width x depht (dyna)
311156-99-80-0	311156-99-60-0	311156-99-90-0	430x360	420x370
306735-99-80-0	306735-99-60-0	306735-99-90-0	430x410	450x420
306431-99-80-0	306431-99-60-0	306431-99-90-0	430x460	450x500
311398-99-80-0	311398-99-60-0	311398-99-90-0	430x510	450x530
306471-99-80-0	306471-99-60-0	306471-99-90-0	430x560	450x560
306470-99-80-0	306470-99-60-0	306470-99-90-0	480x460	500x500
311401-99-80-0	311401-99-60-0	311401-99-90-0	480x510	500x530
306472-99-80-0	306472-99-60-0	306472-99-90-0	480x560	500x560

COMFORT CUSHIONS CORPUS

Check	Mesh	Width x depht (sits)	Width x depht (dyna)
101499-99-0	104739-99-0	430x410	450x420
101225-99-0	104725-99-0	430x460	450x500
101511-99-0	104741-99-0	430x510	450x530
101516-99-0	104742-99-0	430x560	450x560
101517-99-0	104743-99-0	480x460	500x500
101519-99-0	104744-99-0	480x510	500x530
101518-99-0	104745-99-0	480x560	500x560

VS seat

A seat with a standing function that can be both useful and practical. ICS (Intelligent Control System) gives you single button raising from all positions and the option to save your favorite positions. You are raised smoothly and stand steadily on the footplate, which descends to the ground

during raising. You can choose a seat that you can drive standing up. The footplate then has slightly larger wheels. There is a wide range of accessories for the seat, including the UniTrack accessory system.



Arm rest bag
100671-99-0 Right
100672-99-0 Left
102588-99-0 Front of seat frame



Parallel movement panel holder
104302-99-0



Table, adjustable angle
100575-99-0 (350 mm)
100576-99-0 (400 mm)
100577-99-0 (450 mm)



Small table, incl. bracket
100578-99-0 Fitted on the right
101691-99-0 Fitted on the left



Table, swiveling with central cutout
100573-99-0 (600 mm)
100574-99-0 (650 mm)

With R-net Electronics
103980-99-0
103981-99-0



Bag hook set
101473-99-0



Belt buckle, manual adjustment
100518-99-0



Table, angle adjustable
103426-99-0



Table, swiveling with side cutout
100571-99-0 (600 mm)
100572-99-0 (650 mm)

With R-net Electronics
103982-99-0
103983-99-0




Mirror, incl. bracket
101440-99-0



Co-Pilot
103184-99-0

With R-net Electronics
104002-99-0

VS seat



Back rest reinforcement
101101-99-0



Chest roll, adjustable
101099-99-0 Imitation leather,
340 mm back
101095-99-0 Imitation leather,
390 mm back
101097-99-0 Imitation leather,
440 mm back



Foot strap with plastic lock
103427-99-0



Side support
101261-99-0 Right outside
101262-99-0 Left outside



Adjustable head rest
103460-99-0

Lateral headrest adjustment
103697-99-0



Arm rest raising set
100523-99-0



Calf rest, Sr.
103394-99-0

Calf rest, Jr.
104137-99-0



Push handle set
101100-99-0



Manual footplate fold/lock
100979-99-0



Foot strap with padding
103104-00-0



Side support
101263-99-0 Right inside
101264-99-0 Left inside



Head rest, adjustable
104130-99-0 Fabric
104131-99-0 Imitation leather

Lateral head rest adjustment
103697-99-0



Inflatable support for the curve of the back
101972-99-0



Mug/Phone holder
103044-99-0



Sliding backrest
104115-99-0 (340x620 mm)
104116-99-0 (390x620 mm)
104117-99-0 (440x620 mm)

The calf rest is specific to the VS seat.
For other supports, see the UniTrack
accessory system.

Conversion articles, VS seat



Arm Rest Neutral, Fabric

1822782, 250 mm
1822783, 300 mm
1822784, 350 mm



Arm Rest Neutral, Imitation Leather

1822786, 250 mm
1822787, 300 mm
1822788, 350 mm



Chest roll, imitation leather

306110-99-90-0 For 340 mm back
306108-99-90-0 For 390 mm back
302956-99-90-0 For 440 mm back



Shaft, arm rest bracket

306482-99-0 For 340 mm back
302859-99-0 For 390 mm back
304969-99-0 For 440 mm back



Back set, 550 mm, high

310291-00-0 (340 mm wide)
310251-00-0 (390 mm wide)
310252-00-0 (440 mm wide)



Back set, 620 mm, high

305480-00-UD-0 (340 mm wide)
306515-00-UD-0 (390 mm wide)
306490-00-UD-0 (440 mm wide)



Seat widening

1821779 (500 mm) for VS Jr
1821774 (550 mm) for VS Sr



Seat extension set

303006-28-UD-0

Seat and back cushions, VS

The dimensions below are for the seat. The actual dimensions of the cushion are normally a little bigger.



BACK CUSHIONS VS

Mesh	Imitation leather	width x height (seat)	width x height (cushion)
307506-99-60-0	307506-99-90-0	340x550	370x530
307507-99-60-0	307507-99-90-0	390x550	420x530
307508-99-60-0	307508-99-90-0	440x550	470x530
306475-99-60-0	306475-99-90-0	340x620	370x600
306430-99-60-0	306430-99-90-0	390x620	420x600
306473-99-60-0	306473-99-90-0	440x620	470x600

SEAT CUSHIONS VS

Mesh	Imitation leather	width x depht (seat)	width x depht (cushion)
311156-99-60-0	311156-99-90-0	430x360-400 1,2 A	420x370
306735-99-60-0	306735-99-90-0	430x425 B	450x420
311397-99-60-0	311397-99-90-0	430x450 C	450x450
306431-99-60-0	306431-99-90-0	430x475-525 D, E, F	450x500
311398-99-60-0	311398-99-90-0	430x510	450x530
306471-99-60-0	306471-99-90-0	430x550 G	450x560
311399-99-60-0	311399-99-90-0	480x450 C	500x450
306470-99-60-0	306470-99-90-0	480x475-500 D, E	500x500
311401-99-60-0	311401-99-90-0	480x510	500x530
306472-99-60-0	306472-99-90-0	480x525-550 F, G	500x560

COMFORT CUSHIONS VS

Check	Mesh	width x depht (seat)	width x depht (cushion)
101499-99-0	104739-99-0	430x410 B	450x420
101508-99-0	104740-99-0	430x460 C	450x450
101225-99-0	104725-99-0	430x475-500 D, E	450x500
101511-99-0	104741-99-0	430x525 F	450x530
101516-99-0	104742-99-0	430x550 G	450x560
101517-99-0	104743-99-0	480x475-500 D, E	500x500
101519-99-0	104744-99-0	480x525 F	500x530
101518-99-0	104745-99-0	480x550 G	500x560

The letter shown in the cushion name defines wich Seat Depth the cusion fits.

More information about adjusting seas sizes can be found in the seating system's Service Manual.

PS seat

A seat that is very reminiscent of manual wheelchair seats. You adjust the seat with simple tools. With a sliding, telescopic frame and expanding seat pan, you can change the size


and height when you want. The level of adjustability is unmatched. The armrests both flip up and move out of the way and the footplate hangers can swing away or swing in.



Belt buckle, manual adjustment
100518-99-0



Belt, self-adjusting
100519-99-0 Standard
100600-99-0 Long



Bag hook set
101473-99-0



Co-Pilot
103184-99-0
With R-net Electronics
104001-99-0



Adjustable head rest (adjustable back)
100903-99-0



Trunk support
103056-99-0



Head rest, adjustable (Corpus back)
104130-99-0 Fabric
104131-99-0 Imitation leather



Table angle adjustable
103426-99-0 (500 x 332 mm)



Back box
102059-99-0



Mug/Phone holder
103044-99-0



Bracket UniTrack
103680-99-0

TS seat

Robust, flexible seat in a bold design with neat solutions. The backrest and seat width can both be adjusted as required. Ergonomic with sturdy Velcro fastenings on the back and seat surfaces, the whole seat can be moulded to body contours. The seat frame has an integrated UniTrack rail to allow stepless adjustment of the seat depth and armrest height.



Seat pan
104736-99-0



Parallel movement armrest kit (for powered recline)
(med elrygg)
104348-99-0



Parallel movement panel holder
104302-99-0



Crutch holder
101956-99-0



Oxygen holder
101964-99-0



Calf rest
102962-99-0



Belt buckle, manual adjustment
100518-99-0



Belt self adjusting
100519-99-0



Belt self adjusting, long
100600-99-0



Push handle
101956-99-0



Headrest
104361-99-0



Trunk supports
103056-99-0



Thigh supports
103095-99-0

Conversion articles, TS



Backrest height extension kit
(150 mm)
104348-99-0



Seat cushions
1823716, Smal
1823717, Bred



Tension straps (seat)
104348-99-0



Backrest cushions
1823714, Smal
1823715, Bred



Tension straps (backrest)
1823725, Låg rygg
1823723, Hög rygg

CS seat

A seat for those who can cope with simple functions. It has a slightly lower back rest and manually adjustable back rest angle but the comfort required to allow you to sit comfortably for an extended period of time.



Belt buckle, manual adjustment
100518-99-0



Belt, self-adjusting
100519-99-0 Standard
100600-99-0 Long



Head rest
102149-99-0



Mug/Phone holder
103044-99-0



Crutch holder
101956-99-0



Oxygen holder
101964-99-0



Foot strap with padding
103104-00-0



Foot strap with plastic lock
103427-99-0

Miniflex seat

A child seat for children up to 35 kg. It is easy to use so that the child immediately feels safe. For smaller children, it fits best on one of our child chassis. For larger children, on the

other hand, it is an advantage to combine the seat with an adult chassis. This makes it easy to switch to an adult seat when the time comes.



Head rest
100748-99-0 Fabric
100749-99-0 Imitation leather



Side support
100750-99-0 Fabric
100751-99-0 Imitation leather



Thigh rest
100752-99-0 Fabric
100753-99-0 Imitation leather



Knee stop
100866-99-0 Fabric
100867-99-0 Imitation leather



Belt buckle, manual adjustment
100518-99-0



Four-point belt, incl. bracket
100598-99-0



Parallel movement panel holder
101046-99-0 Right
101047-99-0 Left



Belt, self-adjusting
100519-99-0



Foot strap with plastic lock
103427-99-0



Foot strap with padding
103104-00-0



Mug/Phone holder
103044-99-0



Co-Pilot
103182-99-0

Conversion articles, Miniflex



Arm rest, 200

308021-00-30-1 Right, fabric
308021-00-30-2 Left, fabric



Arm rest, 200

308021-00-90-1 Right, imitation leather
308021-00-90-2 Left, imitation leather



Arm rest, 250

307400-00-30-1 Right, fabric
307400-00-30-2 Left, fabric



Arm rest, 250

307400-00-90-1 Right, imitation leather
307400-00-90-2 Left, imitation leather



Seat plate set

306343-28-30-0 (250x250)
306343-28-90-0 (250x250)
306344-28-30-0 (300x300)



Seat plate set

306344-28-90-0 (300x300)
306345-28-30-0 (350x350)
306345-28-90-0 (350x350)



Seat cushion

306384-99-30-0 (250x250)
306384-99-90-0 (250x250)
306385-99-30-0 (300x300)



Seat cushion

306385-99-90-0 (300x300)
306386-99-30-0 (350x350)
306386-99-90-0 (350x350)



Back set

306351-28-30-0 (250x350)
306351-28-90-0 (250x350)
306351-28-30-0 (250x350)



Back set

306351-28-90-0 (250x350)
306353-28-30-0 (350x500)
306353-28-90-0 (350x500)



Back cushion

306381-99-30-0 (250x350)
306381-99-90-0 (250x350)
306382-99-30-0 (300x400)



Back cushion

306382-99-90-0 (300x400)
306383-99-30-0 (350x500)
306383-99-90-0 (350x500)



footplate

307962-00-0 Small
307963-00-0 Large

Our new accessory system - UniTrack

The system consists of a few simple basic parts that you combine and fix in a track on the chair. Then choose the cushion you want to have and you get support in exactly the right place with the right shape. This creates infinite possibilities for adaptation and adjustment.

Accessory tracks

The track, the basis of our new accessory system, has a strategic location on the seat and the rear of the cushion to provide fixing points. The track allows infinitely variable adjustment. All parts of the system have standard dimensions of 25 mm center-to-center and are fixed with standard 6 mm screws.

Adjustment elements

These parts make it possible to rotate the support's cushions although the fixing point remains fixed. The parts are incredibly strong and adjustable in such small steps that the result is almost infinitely variable. Combine and place the support in exactly the right place.

Quick release fitting

A half turn of the handle is enough.
For the trunk rest it is also possible to fit a swing away lock.

Support cushions

Cushions in several sizes and shapes with covers made of breathable material. The cushions are easy to attach to or remove from other system parts. The cushions are given only single letter designations to avoid a cushion being associated with a particular area of use.

The letter and size will be included in the designation of the support.
Example: if you want to order a calf rest and use cushion A, the designation will be Calf rest A 150x140.



A	150x140		
B	110x90		
C	200x90		
D	200x200		
E	220x150		
F	220x90		
G	300x90		
H	170x175		
J	120x165		

Put them together and you can create the following, among other things:

THIGH REST/HIP REST



Adjustment of depth using tracks, adjustment of height using stop screws.



Adjustment of angle.



Adjustment of angle and width.

CALF REST



Thanks to the channels, the cushion can be rotated 360°.



Adjustment of angle and depth.



Adjustment of height and width.

AMPUTATION SUPPORT



Easy to remove using the quick release fitting.



Adjustment of angle and height.



Adjustment of depth, width and angle.



Can be fitted at a 30-degree angle and the cushion can be shaped.

TRUNK REST



Adjustment of height.



Adjustment of angle and width.



Adjustment of height using stop screws.



Adjustment of depth.

Support, UniTrack



Calf rest B 110x90
103055-99-0



Calf rest J 165x120
102982-99-0



Calf rest A 150x140
103058-99-0



Calf rest H 170x175
102978-99-0



Calf rest D 200x200
103060-99-0



Thigh support B 110x90
103092-99-0



Thigh support J 165x120
103094-99-0



Thigh support C 200x90
103095-99-0



Thigh support G 300x90
102977-99-0



Trunk support A 150x140
103050-99-0



Trunk support F 220x90
103056-99-0



Trunk support E 220x140
103062-99-0

Support, UniTrack



Hip rest B 110x90
103077-99-0



Hip rest A 150x140
103079-99-0



Hip rest C 200x90
103081-99-0



Hip rest D 200x200
103083-99-0



**Amputation support A
150x140**
103067-99-0



**Amputation support J
165x120**
103068-99-0



**Amputation support H
170x175**
103070-99-0



**Amputation support D
200x200**
103071-99-0

Spare parts, UniTrack



Amputation support
1821686



Calf rest
1821681



Trunk support, left
1821684



Trunk support, right
1821685



Hip rest
1821682



Thigh support
1821683



Cushion Model A
1821689



Cushion Model B
1821690



Cushion Model C
1821695



Cushion Model D
1821691



Cushion Model E
1821692



Cushion Model F
1821693



Cushion Model G
1821696



Cushion Model H
1821694

Spare parts, UniTrack



Cushion Model J
1821245



Cover Model B
1821698



Cover Model D
1821700



Cover Model F
1821702



Cover Model H
1821704



**Track, Corpus,
430 mm width**
1822380 (460 mm depth)
1822381 (510 mm depth)
1822382 (560 mm depth)



**Track, VS,
430 mm width**
1821772 for VS Jr
1821773 for VS Sr



Cover Model A
1821697



Cover Model C
1821699



Cover Model E
1821701



Cover Model G
1821703



Cover Model J
1821706



**Track, Corpus,
480 mm width**
1822383 (460 mm depth)
1822384 (510 mm depth)
1822385 (560 mm depth)



**Track, VS,
480 mm width**
1821779 for VS Jr
1821774 for VS Sr

Control system P+/R-net/VR2



Table with flip over joystick (Rnet only)

Corpus back 390 mm
103974-99-0

Corpus back 440 mm
103975-99-0



Table with flip over joystick (Rnet only)

VS back 390 mm
103978-99-0

VS back 440 mm
103979-99-0



Kompakt joystick (Rnet only)

103817-99-0



Kompakt joystick Advanced (Rnet only)

103814-99-0



Chin Control Kompakt joystick (Rnet only)

103817-99-0

Switch set
104641-99-0

Holder
104143-99-0



Attendant Control, (Rnet only)

104222-99-0



Push button switch box R-net

102798-99-0, buttons

102799-99-0, toggle switch



Pc Mouse Control (Rnet only)

104688-99-0



Attendant Control, Pilot+

101205-99-0

Panel holder extra
100958-99-0



Extra switch 3-4 Power functions VR2

101623-99-0



Rocker set, Pilot+

100524-99-0



Panel holder back rest

101204-99-0



Articulated panel holder Central control, Pilot+

101260-99-0

Control system P+/R-net/VR2



Joystick grip, ball
100962-99-0



Joystick grip, pin
100965-99-0



Joystick grip, ball PG8/P+
100963-99-0



Joystick grip, mushroom
100964-99-0



Joystick grip, T
100961-99-0



Joystick grip, funnel
100960-99-0

DriveLink MultiController



The MultiController is the interface that allow you to control your environment from your wheelchair.

Problems are there to be solved. This is the idea behind our alternative control system, DriveLink. The system connects a number of functions directly to the basic unit, which makes it possible for you to use various solutions to control your wheelchair.

The Multi Controller has a number of built in functions as standard, for example control of home entertainment appliances (TV, Video, DVD) via your chosen input device. Accessories for the Multi Controller include, among other, things a mobile phone module, but it also acts as an interface to other aftermarket environmental controls, allowing you to activate computers, IR door openers and other functions, from your wheelchair.

DriveLink MC

- Choice of control system
- Seat functions (max. 9)
- Customized display
- Right and left mouse click
- Functions independently
- Recordable message system
- Two active control systems
- 18 possible outputs
- IR mouse and IR games
- Environment control
- Voice-controlled feedback
- Telephony





Drive Link Multicontroller Basic Control Unit

I04615-99-0, Pilot+
I03822-99-0, R-net

Includes:

Electronics	612899-99-0, Pilot + 1822779, R-net
Display	612900-00-0
Cable	309699-00-0
Holder	1821583
Bag	1821541

Please note: Complete Set= 1 x Basic Control Unit + 1 Input Device

INPUT DEVICES

Joystick control
I03823-99-0

Joystick
Holder display



Mini-joystick
I03824-99-0



Chin Control
Kompakt joystick
I03895-99-0

Switch set
I04641-99-0
Holder
I04143-99-0



Chin control
I03827-99-0



Finger control
I03826-99-0



Head control
I03828-99-0

3-5 Switch head control
Head support
Bracket head support



Sip and puff
I03829-99-0

Siff and puff unit
Stand, chin control/SoBV.2



Tablet Control
I04128-99-0



Multi button control
I03830-99-0

Contact box



DriveLink Accessories



Holder, Sip and puff/Chin control el.
101210-99-0



DL IR-receiver
101215-99-0



DL+ Gyro unit
101212-99-0



Programming kit for Multicontroller
102568-99-0



DL Multimedia for Multicontroller
102561-99-0



DL Easy Mobile
102555-99-0



DL+ Function switch touch
101207-99-0



DL Function switch
101222-99-0



Holder Chin control
104143-99-0



External on/off
101208-99-0




Table for Tablet Control
1821941





Holder Mini-joystick
104144-99-0




DL Holder GEWA
101216-99-0

Permolock C

Safety and freedom are two of our key words when we develop our wheelchair solutions. An electric wheelchair increases your mobility in itself, but for real freedom it must also be safe and easy to transport the chair in a vehicle. All our chairs have mounting eyes so that you can secure them with straps during transport. For our combination chairs, C400 and C500, there is now also a very safe docking lock for vehicles, the Permolock C. The Permolock C has been subjected to extensive safety tests, including crash tests at Volvo's test facility. With the Permolock C you can sit in your wheelchair as safely as in the vehicle's normal seats. The lock can be installed in the driver's seat or on the passenger side and can be switched between them as required.




Permolock C
102840-99-0



Chair base
1820550



Plunge
102842-99-0



Quick anchoring set
1820549



Installation spacer
1820547



Quick release fitting release handle
1820548



Release handle
314036-99-0

Permolock C3


The Permolock C3 is the latest addition to the Permolock family, designed specifically for the C300, C350, K300 and K450 models. Thanks to its compact format, the Permolock C3 is easy to mount inside a car. The C3 is adjustable in height and angle which makes mounting easy, even on sloping car floors. Its smart design makes it service-friendly, requiring only the removal of the top cover to access vital parts. A new function for the C3 is that it allows the user to drive into or, reverse into, the lock, depending on where it is mounted in the car. The lock can be mounted in the driver's position or the passenger positions (front or rear) in cars suitable for transport of powered wheelchairs.




Permolock C3
1823300 with switch std.
1824756 with switch large



Release handle
314036-99-0
Quick release fitting release handle
1820548




Switch std.
1820525




Locking pins, Permolock C3
1824229




Switch large
1824735



Chair base
1820545
Extender chair base
1821551




C350 Reinforcement
1824736



Quick anchoring set
1820549



Installation spacer
1824737



Extension cable
1823604